

B-Lite® MESMO® Mème® RR

Moderate Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
15625-130LW25¹	93	93	32	59	59	130
15625-150LW25¹	98	98	34	62	62	150
15625-175LW25¹	103	103	35	65	65	175
15625-200LW25¹	108	108	36	68	68	200
15625-230LW25¹	113	113	37	71	71	230
15625-260LW25¹	118	118	39	74	74	260
15625-290LW25¹	123	123	40	77	77	290
15625-325LW25¹	128	128	41	80	80	325
15625-365LW25¹	133	133	42	83	83	365
15625-405LW25¹	138	138	44	86	86	405
15625-450LW25¹	143	143	45	89	89	450
15625-495LW25¹	148	148	46	92	92	495
15625-545LW25¹	153	153	47	95	95	545
15625-600LW25¹	158	158	49	98	98	600

Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
15626-165LW25¹	93	93	39	64	64	165
15626-190LW25¹	98	98	41	67	67	190
15626-215LW25¹	103	103	43	71	71	215
15626-250LW25¹	108	108	44	74	74	250
15626-285LW25¹	113	113	46	77	77	285
15626-325LW25¹	118	118	48	80	80	325
15626-365LW25¹	123	123	49	83	83	365
15626-410LW25¹	128	128	51	87	87	410
15626-460LW25¹	133	133	53	90	90	460
15626-510LW25¹	138	138	54	93	93	510
15626-565LW25¹	143	143	56	96	96	565
15626-625LW25¹	148	148	58	100	100	625
15626-690LW25¹	153	153	59	103	103	690
15626-760LW25¹	158	158	61	106	106	760

Extra Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
15627-195LW25¹	93	93	47	69	69	195
15627-225LW25¹	98	98	49	73	73	225
15627-260LW25¹	103	103	51	77	77	260
15627-300LW25¹	108	108	53	80	80	300
15627-345LW25¹	113	113	55	84	84	345
15627-390LW25¹	118	118	57	87	87	390
15627-440LW25¹	123	123	59	91	91	440
15627-495LW25¹	128	128	61	94	94	495
15627-555LW25¹	133	133	63	98	98	555
15627-615LW25¹	138	138	65	101	101	615
15627-685LW25¹	143	143	67	105	105	685
15627-755LW25¹	148	148	69	108	108	755
15627-835LW25¹²	153	153	71	112	112	835
15627-920LW25¹²	158	158	73	115	115	920

B-Lite® MESMO® Replicon® AR

Moderate Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
15635-155LW25¹	103	103	35	50	69	155
15635-175LW25¹	108	108	36	52	72	175
15635-200LW25¹	113	113	37	55	75	200
15635-230LW25¹	118	118	39	57	78	230
15635-255LW25¹	123	123	40	59	82	255
15635-290LW25¹	128	128	41	61	85	290
15635-320LW25¹	133	133	42	63	88	320
15635-355LW25¹	138	138	44	66	91	355
15635-395LW25¹	143	143	45	68	94	395
15635-435LW25¹	148	148	46	70	97	435
15635-480LW25¹	153	153	47	72	100	480
15635-525LW25¹	158	158	49	74	103	525

Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
15636-190LW25¹	103	103	43	57	76	190
15636-220LW25¹	108	108	44	59	79	220
15636-250LW25¹	113	113	46	62	83	250
15636-285LW25¹	118	118	48	64	86	285
15636-320LW25¹	123	123	49	67	89	320
15636-360LW25¹	128	128	51	69	93	360
15636-405LW25¹	133	133	53	72	96	405
15636-450LW25¹	138	138	54	74	100	450
15636-500LW25¹	143	143	56	77	103	500
15636-550LW25¹	148	148	58	79	107	550
15636-610LW25¹	153	153	59	82	110	610
15636-670LW25¹	158	158	61	84	114	670

Extra Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
15637-200LW25¹	98	98	49	61	79	200
15637-230LW25¹	103	103	51	63	83	230
15637-265LW25¹	108	108	53	66	87	265
15637-300LW25¹	113	113	55	69	90	300
15637-345LW25¹	118	118	57	72	94	345
15637-385LW25¹	123	123	59	75	98	385
15637-435LW25¹	128	128	61	78	102	435
15637-485LW25¹	133	133	63	80	105	485
15637-545LW25¹	138	138	65	83	109	545
15637-605LW25¹	143	143	67	86	113	605
15637-665LW25¹	148	148	69	89	117	665
15637-735LW25¹	153	153	71	92	121	735
15637-810LW25¹²	158	158	73	95	124	810

B-Lite® MESMO® Opticon® AS

Moderate Projektion


Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
15645-170LW25¹	113	95	37	51	67	170
15645-190LW25¹	118	100	39	53	70	190
15645-215LW25¹	123	104	40	55	72	215
15645-240LW25¹	128	108	41	57	75	240
15645-270LW25¹	133	112	42	59	78	270
15645-300LW25¹	138	116	44	61	81	300
15645-330LW25¹	143	121	45	63	83	330
15645-365LW25¹	148	125	46	65	86	365

Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
15646-210LW25¹	113	95	46	58	75	210
15646-240LW25¹	118	100	48	61	78	240
15646-270LW25¹	123	104	49	63	81	270
15646-305LW25¹	128	108	51	65	84	305
15646-340LW25¹	133	112	53	68	87	340
15646-380LW25¹	138	116	54	70	90	380
15646-420LW25¹	143	121	56	72	93	420
15646-465LW25¹	148	125	58	75	96	465
15646-510LW25¹	153	129	59	77	99	510

Extra Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
15647-290LW25¹	118	100	57	69	86	290
15647-325LW25¹	123	104	59	71	89	325
15647-365LW25¹	128	108	61	74	93	365
15647-410LW25¹	133	112	63	77	96	410



B-Lite® Microthane® Même® RR

Moderate Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
30625-130LW25 ¹	96	96	34	62	62	130
30625-175LW25 ¹	106	106	37	68	68	175
30625-200LW25 ¹	111	111	38	71	71	200
30625-230LW25 ¹	116	116	39	74	74	230
30625-260LW25 ¹	121	121	41	77	77	260
30625-290LW25	126	126	42	80	80	290
30625-325LW25	131	131	43	83	83	325
30625-365LW25	136	136	44	86	86	365
30625-405LW25	141	141	46	89	89	405
30625-450LW25 ¹	146	146	47	92	92	450

Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
30626-190LW25 ¹	101	101	43	71	71	190
30626-215LW25 ¹	106	106	45	74	74	215
30626-250LW25	111	111	46	77	77	250
30626-285LW25	116	116	48	80	80	285
30626-325LW25	121	121	50	83	83	325
30626-365LW25	126	126	51	87	87	365
30626-410LW25	131	131	53	90	90	410
30626-460LW25	136	136	55	93	93	460
30626-510LW25	141	141	56	96	96	510
30626-565LW25 ¹	146	146	58	100	100	565
30626-625LW25 ¹	151	151	60	103	103	625
30626-690LW25 ¹	156	156	61	106	106	690

Extra Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
30627-260LW25 ¹	106	106	53	80	80	260
30627-300LW25	111	111	55	84	84	300
30627-345LW25	116	116	57	87	87	345
30627-390LW25	121	121	59	91	91	390
30627-440LW25	126	126	61	94	94	440
30627-495LW25	131	131	63	98	98	495
30627-555LW25	136	136	65	101	101	555
30627-615LW25	141	141	67	105	105	615
30627-685LW25	146	146	69	108	108	685
30627-755LW25 ¹	151	151	71	112	112	755
30627-835LW25 ^{1,2}	156	156	73	115	115	835

B-Lite® Microthane® Replicon® AR

Moderate Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
30635-155LW25 ¹	106	106	37	52	72	155
30635-175LW25 ¹	111	111	38	55	75	175
30635-200LW25 ¹	116	116	39	57	78	200
30635-230LW25 ¹	121	121	41	59	82	230
30635-255LW25 ¹	126	126	42	61	85	255
30635-290LW25 ¹	131	131	43	63	88	290
30635-320LW25 ¹	136	136	44	66	91	320
30635-355LW25 ¹	141	141	46	68	94	355
30635-395LW25 ¹	146	146	47	70	97	395
30635-435LW25 ¹	151	151	48	72	100	435
30635-480LW25 ¹	156	156	49	74	103	480

Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
30636-190LW25 ¹	106	106	45	59	79	190
30636-220LW25 ¹	111	111	46	62	83	220
30636-250LW25 ¹	116	116	48	64	86	250
30636-285LW25	121	121	50	67	89	285
30636-320LW25	126	126	51	69	93	320
30636-360LW25	131	131	53	72	96	360
30636-405LW25	136	136	55	74	100	405
30636-450LW25	141	141	56	77	103	450
30636-500LW25	146	146	58	79	107	500
30636-550LW25 ¹	151	151	60	82	110	550
30636-610LW25 ¹	156	156	61	84	114	610
30636-670LW25 ¹	161	161	63	87	117	670

Extra Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
30637-230LW25 ¹	106	106	53	66	87	230
30637-265LW25 ¹	111	111	55	69	90	265
30637-300LW25 ¹	116	116	57	72	94	300
30637-345LW25	121	121	59	75	98	345
30637-385LW25	126	126	61	78	102	385
30637-435LW25	131	131	63	80	105	435
30637-485LW25	136	136	65	83	109	485
30637-545LW25 ¹	141	141	67	86	113	545
30637-605LW25 ¹	146	146	69	89	117	605
30637-665LW25 ¹	151	151	71	92	121	665
30637-735LW25 ¹	156	156	73	95	124	735
30637-810LW25 ^{1,2}	161	161	75	97	128	810

B-Lite® Microthane® Opticon® AS

Moderate Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
30645-170LW25 ¹	116	98	39	53	70	170
30645-190LW25 ¹	121	103	41	55	72	190
30645-215LW25 ¹	126	107	42	57	75	215
30645-240LW25 ¹	131	111	43	59	78	240
30645-270LW25 ¹	136	115	44	61	81	270
30645-300LW25 ¹	141	119	46	63	83	300
30645-330LW25 ¹	146	124	47	65	86	330
30645-365LW25 ¹	151	128	48	67	89	365

Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
30646-210LW25 ¹	116	98	48	61	78	210
30646-240LW25 ¹	121	103	50	63	81	240
30646-270LW25	126	107	51	65	84	270
30646-305LW25	131	111	53	68	87	305
30646-340LW25	136	115	55	70	90	340
30646-380LW25	141	119	56	72	93	380
30646-420LW25 ¹	146	124	58	75	96	420
30646-465LW25 ¹	151	128	60	77	99	465
30646-510LW25 ¹	156	132	61	80	102	510

Extra Hohe Projektion

Bestellnr.	Dimensionen (mm)					Vol. (ml)
	A	B	C	D	LVC	
30647-290LW25 ¹	121	103	59	71	89	290
30647-325LW25 ¹	126	107	61	74	93	325
30647-365LW25 ¹	131	111	63	77	96	365
30647-410LW25 ¹	136	115	65	79	100	410