



Niedrige Projektion					SLS	
Bestellnr.	Dimensionen (mm)				LVC (mm)	Vol. (ml)
	A	B	C	D		
10724-110 ^{1*}	90	90	26	54	54	110
10724-125 ^{1*}	95	95	27	57	57	125
10724-150 ^{1*}	100	100	28	60	60	150
10724-165 ^{1*}	105	105	29	62	62	165
10724-185*	110	110	29	65	65	185
10724-215*	115	115	30	68	68	215
10724-235*	120	120	31	71	71	235
10724-270 ^{1*}	125	125	32	73	73	270
10724-295 ^{1*}	130	130	33	76	76	295
10724-330 ^{1*}	135	135	34	79	79	330
10724-365 ^{1*}	140	140	35	82	82	365

Moderate Projektion					SMS	
Bestellnr.	Dimensionen (mm)				LVC (mm)	Vol. (ml)
	A	B	C	D		
10725-120 ^{1*}	85	85	32	56	56	120
10725-145 ^{1*}	90	90	33	59	59	145
10725-165 ^{1*}	95	95	35	62	62	165
10725-195*	100	100	36	65	65	195
10725-220*	105	105	37	68	68	220
10725-255*	110	110	38	71	71	255
10725-285*	115	115	40	74	74	285
10725-300*	118	118	40	75	75	300
10725-320*	120	120	41	77	77	320
10725-360*	125	125	43	80	80	360
10725-400*	130	130	44	83	83	400
10725-445*	135	135	45	86	86	445
10725-495*	140	140	47	89	89	495

Hohe Projektion					SHS	
Bestellnr.	Dimensionen (mm)				LVC (mm)	Vol. (ml)
	A	B	C	D		
10726-125 ^{1*}	80	80	37	58	58	125
10726-155 ^{1*}	85	85	39	61	61	155
10726-180 ^{1*}	90	90	40	64	64	180
10726-210 ^{1*}	95	95	42	67	67	210
10726-235*	100	100	44	71	71	235
10726-255*	103	103	45	72	72	255
10726-275*	105	105	45	74	74	275
10726-295*	108	108	46	75	75	295
10726-315*	110	110	47	77	77	315
10726-335*	113	113	48	79	79	335
10726-360*	115	115	49	80	80	360
10726-400*	120	120	51	83	83	400
10726-450*	125	125	53	87	87	450
10726-505*	130	130	55	90	90	505
10726-560*	135	135	57	93	93	560
10726-620 ^{1*}	140	140	59	96	96	620
10726-690 ¹	145	145	61	100	100	690
10726-760 ^{1,2*}	150	150	62	103	103	760
10726-835 ^{1,2}	155	155	64	106	106	835

Extra Hohe Projektion					SXS	
Bestellnr.	Dimensionen (mm)				LVC (mm)	Vol. (ml)
	A	B	C	D		
10727-180 ^{1*}	85	85	47	66	66	180
10727-215 ^{1*}	90	90	49	69	69	215
10727-250 ^{1*}	95	95	51	73	73	250
10727-285*	100	100	53	77	77	285
10727-330*	105	105	55	80	80	330
10727-380*	110	110	57	84	84	380
10727-430*	115	115	59	87	87	430
10727-485*	120	120	61	91	91	485
10727-545*	125	125	63	94	94	545
10727-610*	130	130	65	98	98	610
10727-675	135	135	67	101	101	675
10727-755 ^{1,2*}	140	140	69	105	105	755